


YOSHIDA AIR TURBINE HANDPIECE 600SERIES



Powerful tip-air from 4 spray pots.

Outstanding cooling effectiveness by (2x2 Spray) suppress heat when cutting and cause less pain to patient during treatment.

This handpiece assures efficient water spray into the cavity to remove all debris during cutting.

In addition, it keeps an excellent operation sight.

LINE UP

Fiber Optic Type

High torque 635BP-LUX
Super load 640BP-LUX



Non-Optic Type

High torque 635BQP
Super load 640BQP



YOSHIDA AIR TURBINE HANDPIECE

600SERIES

Why do Japanese Dentists Choose our turbines? These are the main reasons.

1. Well-Designed features
2. Effective wide-range spraying
3. Strong cutting efficiency
4. Excellent cooling system
5. Provides excellent operation sight
6. User-friendly
7. High Quality



Quick coupling provides 360° rotation flexibility.
Yoshida turbine is easy to operate like a Pen.



- Superior cooling effect -
2x2 Spray suppresses friction, therefore causes less damage of pulp and less pain of patient during the teeth cutting. And it can also provide effective cutting with appropriate spraying system.

- Excellent Operation Sight for the user -
Powerful spraying can target the certain cavity spot and effectively remove the cutting debris by strong spray power. Thus, it gives the users good sight for their treatments, and it also provides the ideal working site without cutting debris.

With YOSHIDA SPRAY, you can lubricate inside the chuck from the bur attachment, as well.

 Sold separately



ALL PRODUCTS are made in JAPAN!

Contact us:

THE YOSHIDA DENTAL MFG. CO.,LTD.
Headquarter & Tokyo factory

1-3-6 Kotobashi, Sumida-ku Tokyo, Japan ZIP130-8516

TEL: +81-3-3631-2165 (International Business Div.)

FAX: +81-3-3631-2685 (International Business Div.)

Ibaraki factory

3262-4 Kamata, Yoshiwara, Amimachi, Inashiki-Gun, Ibaraki, Japan ZIP300-1155

 **YOSHIDA**