

YOSHIDA
Touching your heart

SEASONS





春の桜、夏の花火、秋の月、冬の雪。四季の美しさを、SEASONSの診療室に再現します。

SEASONS adds luxury to your treatment space in a way prime food in season, richest in taste and texture adds luxury to your dish. In Japan, seasonal food at its best is called "Shu-n" (旬), and it is valued as extreme luxury. That's the idea behind SEASONS. SEASONS offers best quality while giving you the freedom to customize your dental unit so you get to have the best work station precisely suited to your needs.

Classic Line

With SEASONS, less is more. You get only the features you need, so your dental equipment is perfectly suited to the treatments you offer. Enjoy Japanese reliability and build quality with SEASONS, while offering warm hospitality to patients through new levels of ergonomic comfort.



SEASONS





SEASONS is packed with patient-friendly features that make dental treatment a comfortable and relaxing experience.

Comfortable patient position



- 1** Nose-knee-leveled posture
Keeping the nose and the knee level while lowering the belly keeps the patient relaxed and comfortable.
- 2** Patient safety
A seat-height that adjusts as low as 380mm – caters for all types of patient, especially small children and elders. The contoured 'bucket' chair also cushions the patients and keeps them in place.

Easy-to-use doctor's table



- 3** Gently curved corners protect users' hands and provide smooth access to instruments.
- *Instruments in the photo are not include in the standard feature.

Easy-to-maintain assistant element




- 4** A swivel feature enables flexible positioning for easy access by assistant.
- *Instruments and panel switch in the photo are not included in the standard feature.



Need more features and a feeling of supreme comfort? Our SEASONS in Luxury Line makes patients feel almost weightless.

SEASONS
旬
OPTIONS



有田
ARITA
Porcelain

匠
TAKUMI
Seat

- ▼ Upgradable factory option
- 1 DOCTOR'S TABLE
- 2 HEADREST
- 3 ARMREST
- 4 UNIT
- 5 ASSISTANT ELEMENT
- 6 FOOT CONTROLLER
- 7 OPERATING LIGHT

▼ Recommended option

*For more details, see the 'Recommended Options' sheet.



Make patients feel like VIPS.


Your chair and spittoon bowl can be upgraded to add more luxury, you can choose your own colors and designs to put patients at their ease.

TAKUMI seat for total relaxation



L-2536


3-layer structure




Luxury-soft sponge for deep comfort

Soft sponge spreads body pressure evenly

Moderately firm sponge supports weight



[Conventional seat]
Body-pressure distribution



[TAKUMI seat]
Body pressure is widely spread



墨色
Rich Black



砥粉色
Fresh Beige


13 TAKUMI seat optional colors

	L-2541		L-8147		L-2540
	L-2543		L-2532		L-2533
	L-8146		L-2536		L-2538
			L-2537		L-2539
			L-8148		L-8149

ARITA Porcelain spittoon bowl for real luxury



Upgrade to our special ARITA porcelain spittoon and enjoy a touch of luxury. Photo-catalytic ceramic coating works as an antibacterial deodorant, too.



Choose your own colors


The Classic Line comes in a variety of upholstery and spittoon colors. Choose the colors that suit your treatment room.



Upholstery




Standard colors




KAGE Black




SAKURA Pink




YUKI White



AI Blue



SORA Blue



WAKABA Green



ICHO Yellow



Optional colors



Metallic Blue



Metallic Bordeaux



Metallic Green



Metallic Black

Spittoon bowl



Standard colors



Clear



Green



Blue



Optional colors



SERiO Clear



SERiO Blue

Upgrade your SEASONS

You can tailor your SEASONS, with some great, upgradeable factory options.



1 DOCTOR'S TABLE



- LED Panel (Bowl rinsing, Lighting, MM Reverse, Optional light function)
- VFD Panel (VFD Display Monitor)
- Anti-Infection System
- Optical circuit
- Instrument circuit
- Table handle (①center ②upper left ③upper right)
- Drive air monometer
- Waste receptacles holder
- Built-in USB receptacle (wiring through the arm)

2 HEADREST



Angled to ensure patient's mouth remains open

Double-articulating head rest

3 ARMREST



4 UNIT



Cup filler and bowl rinsing sensor




Standalone supply box
Allows flexible layout of dental unit

5 ASSISTANT ELEMENT



- Hanger Switch
- Emergency stop switch
- Panel Switch



Assistant hanger is easily moved to side of headrest

6 FOOT CONTROLLER



Chair controller



Step-in type



Sliding Foot Controller




Multi functional foot controller

7 OPERATING LIGHT



Sensor light


Dimensions



Standard equipment



Chair




Gliding headrest



Spittoon bowl



Unit



Operating light (Switch type)



Built-in supply box



Foot controller



Doctor's table*



Assistant element

*Instruments in the photo are not included in the standard.

The product specifications vary depending on the area of purchase. Please contact our international business division for more information.

© THE YOSHIDA DENTAL MFG. CO., LTD.

Distributor

YOSHIDA

Manufacturing site
YOSHIDA SEIKO., LTD.
Address 660 Yakimaki, Namegata-shi, Ibaraki, Japan ZIP 311-3506

THE YOSHIDA DENTAL MFG. CO., LTD.
Address 1-3-6 Kotobashi, Sumida-ku Tokyo, Japan ZIP 130-8516



www.facebook.com/
YoshidaDentalManufacturing



www.yoshida-net.co.jp/en

Second Edition

Recommended Options


'Made in Japan' counts for a lot these days. SEASONS is reliable, durable and exceptionally well built. And now you can customize your dental unit and chair with the options to create your best work station.



DOCTOR'S TABLE



1 Extra stainless tray + mounting hook




Extra table can be installed for more work space (300×258mm)



Extra table for AEH can also be used on the doctor's table (150×190mm)

2 OBON Resin tray (In frost green, frost pink or frost white)



Doctor side
304×408mm (inside)

AEH side
136×176mm (inside)

『What's "OBON"』

Hospitality and OBON tray

OBON is a Japanese tray for serving drinks and for showing hospitality to guests

- 3 Mouse tray
- 4 Film viewer
- 5 Silicon sheet for 1&2

YOSHIDA INSTRUMENTS

※ WITH OPTICAL LIGHT



Royal road SDC1



Super load 640BP-LUX



High torque 635BP-LUX



Super load 640BQP



High torque 635BQP

AIR TURBINE



S40 MI



S40 SI

MICRO MOTOR HANDPIECE



40C

MICRO MOTOR



LX503(LED)



NO.3

QUICK JOINT

DISPLAY ARM



Display arm for monitor
*Monitor in the photo is not included.

UNIT



Water bottle